

## I. Introduction

The use of cell-free DNA (cfDNA) has increasingly become the standard of care for prenatal screening of fetal aneuploidy. NIPT has the potential to yield a non-reportable result. As such, it is important to understand the factors impacting the ability to obtain a clinical result. Maternal weight has an inverse relationship to fetal fraction, while fetal fraction has a direct relationship with gestational age. This study reviews the success rate of obtaining a NIPT result as a function of maternal weight and gestational age (GA).

## II. Methods

A retrospective analysis of 139,995 consecutive maternal blood samples that were submitted to Sequenom Laboratories® (a wholly owned subsidiary of Laboratory Corporation of America® Holdings) for MaterniT® 21 PLUS laboratory developed testing were stratified by maternal weight and gestational age. The percent of NIPT samples that yielded a non-reportable test result was evaluated assuming both factors are independent. Samples were subjected to DNA extraction, library preparation, and whole genome massively parallel sequencing as described by Jensen et al.<sup>1</sup>

## III. Results

Patients <125 lbs have 100% success rate at GA >22 weeks, the lowest success rate is in the highest weight population >300 lbs at 93.7% (ranging from 90.9-100% depending on GA). Stratifying these populations by GA shows only a minor impact in success rate across GA in the <200 lbs population, but a marked improvement as GA increases, matching the average population success rate in the heaviest population at 25 weeks GA.

**Table 1.**  
Success rate by gestational age and maternal weight

Gestational Age	Maternal Weight (lbs)										Average
	<100	100 - 124	125 - 149	150 - 174	175 - 199	200 - 224	225 - 249	250 - 274	275 - 299	>300	
<12 weeks	99.72%	99.71%	99.54%	99.35%	98.97%	97.97%	97.00%	96.11%	95.02%	93.15%	99.03%
13-15 weeks	100.00%	99.85%	99.60%	99.29%	98.98%	98.48%	96.48%	97.65%	94.87%	90.86%	98.94%
16-18 weeks	100.00%	99.70%	99.48%	99.44%	99.14%	98.82%	97.59%	96.99%	95.79%	94.89%	99.06%
19-21 weeks	100.00%	99.62%	99.71%	99.65%	99.38%	98.90%	99.30%	99.02%	99.35%	96.00%	99.45%
22-24 weeks	100.00%	100.00%	99.82%	99.72%	99.73%	99.11%	99.23%	100.00%	100.00%	100.00%	99.70%
25-27 weeks	100.00%	100.00%	99.78%	99.44%	100.00%	99.62%	100.00%	97.40%	100.00%	100.00%	99.66%
28-30 weeks	100.00%	100.00%	99.64%	99.15%	100.00%	100.00%	100.00%	100.00%	100.00%	100.00%	99.70%
>30 weeks	100.00%	100.00%	100.00%	100.00%	100.00%	100.00%	100.00%	100.00%	96.88%	100.00%	99.92%
<b>Average</b>	99.82%	99.74%	99.57%	99.39%	99.09%	98.37%	97.36%	97.05%	95.81%	93.67%	99.09%


**Table 2.**  
Success rate by maternal weight

Maternal Weight (lbs)	Total #	Success Rate
<100	1,339	99.85%
100-124	19,869	99.76%
125-149	42,231	99.57%
150-174	31,941	99.37%
175-199	19,821	99.08%
200-224	11,353	98.29%
225-249	6,706	97.36%
250-274	3,404	96.86%
275-299	1,701	95.71%
>300	1,580	93.48%
<b>Total</b>	139,945	99.09%

**Table 3.**  
Success rate by gestational age

Gestational Age	Total #	Success Rate
<12 weeks	83,416	99.00%
13-15 weeks	25,860	98.90%
16-18 weeks	11,279	99.10%
19-21 weeks	10,276	99.40%
22-24 weeks	4,317	99.70%
25-27 weeks	2,044	99.70%
28-30 weeks	1,335	99.70%
>30 weeks	1,418	99.90%
<b>Total</b>	139,945	99.10%

**Figure 1.** Weekly increase in success rate by maternal weight



## IV. Conclusions

Of the two factors studied, GA and maternal weight, the later has a larger impact on NIPT success rate but it can be improved with an increase in GA. Despite a reduced success rate at extreme maternal weights, especially at early gestational age, cfDNA testing delivers results for more than 93.6% of patients in the >300 lbs population. In this study we show that the non-reportable rate of maternal weight on NIPT results can be improved by waiting to test at a later GA for patients >200 lbs. NIPT can be considered a viable option for aneuploidy screening in obese patients.

## V. References

- Jensen TJ, Zwiefelhofer T, Tim RC, et al. High-throughput massively parallel sequencing for fetal aneuploidy detection from maternal plasma. *PLoS One*. 2013; 8(3):e57381. doi:10.1371/journal.pone.0057381. Epub 2013 Mar 6.